floorline LTC











Our Floorline LTC is the optimal bed solution for home care use or long-term care facilities. Its one piece design makes it easy to fold, store and transport. The Floorline LTC comes in two convenient options including one with a wood-grain finish that compliments any room.

Clinical and Patient Benefits

The Floorline LTC lowers to a distance of 9.5cm/3.7" from the floor. This distance, although seemingly small, has enormous benefits to both the patient and the facility. By just adding 5cm to the distance between the bed and the floor, the impact force of a fall is increased by 50%. This can be the difference between a simple bruise or a fractured hip for the patient that may result in an extented hospital stay or even death. In addition to the physical effects, fall injuries also have psychological effects on the patient as they experience a loss of confidence and fear after a fall.

normally found only on hospital beds, such as Trendelenburg, reverse Trendelenburg and cardiac positioning. With a full assortment of accessories, an easy-to-use handset, and durable LINAK electronics, the Floorline LTC is the best choice for patient care.

Quality

Our products are designed and manufactured to meet all the requirements of relevant standards. At every stage of the development process, every single product is given extensive attention in design, development and manufacturing. This means that a user of this product feels confident with its reliability and the caregiver feels confident with its safety.

Optimal Features

The Floorline LTC features standard positions

HUMAN CARE SWEDEN (HQ) Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com HUMAN CARE AUSTRALIA Unit 4, 49 Bertie Street Port Melbourne, VIC 3207 Phone : +61 1300 132 736 Fax : +61 (03) 9645 4954 info.au@humancaregroup.com HUMAN CARE CANADA 10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613.723.6734 Fax: +613.723.1058 info.ca@humancaregroup.com HUMAN CARE NEDERLAND Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170

info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com

Accessories

ന

IV Pole

Adjustable height. Can be mounted on either end of bed. Maximum SWL 7kg/ 15.4lb. Part no. IVP



Care Assist Rail Designed to assist patients safely in and

out of a bed. Ergonomic design. Part no. 90702



Extension Kit Extends bed by 17.5cm/7" for taller patients. Part no. 90709



Floorline LTC

Handset Five function - included with bed. Part no. 75746 Handset holder Part no. LC2250, Lock Out Key Part no. 75790

Technical Specifications

Self Help Pole Assists patient moving within confines of bed. Maximum SWL 75kg/165.3lb. Part no. SHP600

Side Rails

Available in 1/4 size (Part no. 72790), 1/2 size (Part no. 72791), Split 1/2+1/4 size (Part no. 72792), and 3/4 size (Part no. 72793)

Wall Bumper Bar

Protects the head of the bed. Mounted on lower support bracket between wheels. Part no. 90705



Features

Compact

Folds compactly for storage or transport.



Premium Casters

One central locking and three steering/directional maneuvrability.



Manufactured to all relevant standards including IEC 60601-2-52:2009, EN1970:2002, AS/NZS 3200.2.38:2007

Part no. White Finish	90600 (EU), 90602 (UK), 90604 (NA), 90608 (AU)	Part no. Woodgrain Finish	90601 (EU), 90603 (UK), 90605 (NA), 90609 (AU)
Platform Position (Min-Max)	9.5cm - 70cm / 3.7" - 27.5"	Overall Bed Length (Standard to Extended with Extension Kit)	215cm - 232.5cm / 84.5" - 91.5"
Platform Dimensions (L x W)	215cm x 90.5cm / 84.5" x 35.5"	Platform Dimensions with Extension Kit (L x W)	232.5cm x 90.5cm / 91.5" x 35.5"
Sleeping Dimensions (LxW)	200cm x 90cm / 78.7" x 35.5"	Sleeping Dimensions with Extension Kit (L x W)	217.5cm x 90cm / 85.5" x 35.5"
Output Voltage	Max. 24V DC	IP Rating	IPX4
Input Voltage	120V, 60Hz (90404 to 90607), else 230V, 50Hz	Product Weight	106kg/234lb
Platform Panel Angles	Backrest 70°, Thigh 45°, Knee 110°, Calf 25°, Trendelenburg/ Reverse Trendelenburg 18°	Safe Working Load (SWL)	200kg/440lb
Duty Cycle	2 minutes operation, 18 minutes rest	Maximum User Weight	165kg/363lb
This product is manufactured			

in accordance with all relevant ISO standards.





Follow Us on Social Media





0

in linkedin.com/company/human-care-ab



 $(\mathcal{P}$ pinterest.com/humancaregroup



youtube.com/channel/UC7fyiNh3zTjRx0P1IDf2DaQ

humancaregroup.com

casters for premium

