

Invacare® | **Action® 2NG**



The next generation has arrived

The all-new **Action 2NG** has been developed from the tried-and-tested **Action 3NG** concept.

With new customisation features such as depth-adjustable seat upholstery, the specifications meet the basic needs of the users while offering greater adjustability. The next generation **Action 2NG** shares the same essential platform as the **Action 3NG**, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the **Action 2NG** fits perfectly into the Action family.



Code: BQ152XXXX
(Refer to customer service for the code size you require)

Benefits

Light and comfortable

With a weight of just 15kg, the *Invacare® Action 2NG* fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and the light weight makes transportation trouble-free. Cross-compatibility with other members of the Action family ensures that the *Action 2NG* can be easily adapted to meet individual user needs.

Quality of life

The tried-and-tested technology of the Action family products allow the user to go wherever they desire. Cross-compatibility with *Action 3NG* make this an ideal, basic, fleet product.

Features



NEW Depth-adjustable seat upholstery
Decrease seat depth by 50mm for customised seat size.



Swing in and out footrest
The 80° footrest swings in and out for easy transfers in or out of the chair.



Adjustable center of gravity
The Axle position can be adjusted in height by 5cm.

Technical Data

Warranty:
2 years frame & crossbraces.
Limited 1 year components.

41, 45, 48 cm	40 & 43 cm	80° / 0° - 80°	48,5 cm	7° / 0 - 30°*	43 cm	35 - 48 cm
26 cm	58,5 - 68,5 cm	41 cm	95 - 109 cm	82 / 84,5 cm	93,5 cm	max. 125 kg

* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.



16,2kg Self Propel
14,9kg Transit

Frame Colour



Black

Upholstery



Black Nylon

* Requires reclining back canes.