

# Invacare® | **Action® 3NG**



## **Greater driving performance**

The *Invacare® Action®3NG* has been re-designed to ensure greater driving performance; a lightweight frame and high quality components make for easy rolling and manoeuvrability. An appealing 'next generation' design together with improved functionalities in terms of stability and comfort means that the **Action®3NG** marks the rebirth of our flagship product in style.



**Codes:** CS155XXXX or A3LD-xxXxx L or R  
(Refer to customer service for the code size you require)

## Benefits

### Functionality

Invacare® always seeks to innovate. The **Action®3NG** introduces a new cross bar with a folding mechanism under the seat which improves stability and makes folding and unfolding even easier. Together with a new-style footplate, the Action®3 is a great choice for all-round maximum comfort.

### Durability

Adaptable and modular, the **Action®3NG** accommodates progressive changes in user needs. It is also simple to refurbish, with a wide range of options. This model is a great example of a second-generation wheelchair, with a proven track record behind it.

### Transportable

With a low transfer weight and foldable backrest, the **Action®3NG** is easy to fold and to dismantle, making it easily transportable in a small car.

## Features



### 14kg weight

The total weight of 14kg addresses the need for improved handling, storage and transportation. The lowest transportation weight (without legrest, armrest and rear wheels) is 7.6 kg.



### Brake and Brake extension

Fold the brake handle for easier side transfers.



### Length compensating leg support

New option elevating legrest offers optimal length compensation with the elevation stepped in approx. 3 – 5 degrees.



### Footrest

The proven design of the swing in / swing out footrest has an easy-to-operate quick release mechanism.



### 20 rear wheel positions

Driving performance is enhanced with a choice of 20 rear wheel positions as standard. Users can set the rear wheel to Standard, Medium, Active or Amputee position.



### Magnesium castor attachment

A castor attachment provides ease of angle adjustment and enhanced support. The use of the magnesium attachment reduces the weight of the chair by almost 400 gram.

## Technical Data

### Warranty:

2 years frame & crossbraces.  
Limited 1 year components.



40.5, 45.5, 51 cm



40.5, 45.5, 51 cm



80° / 0° - 80°



38.5, 41, 43.5, 46, 48.5, 51 cm



0° / 7° / 0° - 30°  
/-6° - +24°



40, 43, 51 cm



35 - 48 cm  
34 - 47 cm



20 - 23 - 26 cm,  
24 cm, 19 - 29.5 cm



58.5 - 68.5 cm



sw+ 20.5 cm



100 cm



72 - 85 cm



81.5 - 102 cm



max. 125 kg

\* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.



14.2kg Self Propel  
13.5kg Transit



7.6 kg (without armrest,  
footrest & rearwheels)

### Frame Colour



Anthracite Black

### Upholstery



Black Nylon