

PRESSURE INJURY PREVENTION

Alpha Trancell® Deluxe

Pressure redistribution system



24 hour solution

The Alpha Trancell Deluxe system offers highly effective alternating pressure redistribution in a low cost, simple to use package.

Using a clinically proven 10 minute cycle and a carefully calibrated pressure control on the pump, the Alpha Trancell Deluxe mattress replacement, overlay and seating systems are effective tools for the prevention of pressure ulcers and in the management of superficial tissue damage.

Features

Mattress (Overlay)

- 20 fully alternating modular cells
- Durable polyurethane (PU) material
- Rapid deflate function
- Transport facility
- Securing straps

Mattress (Replacement)

- As above with additional integral foam underlay to provide full depth mattress replacement

Seat Cushion

- 9 fully alternating PU cells
- Modular construction
- Securing straps
- PU coated, knitted fabric cover

Pump

- Manual pressure control with easy to use numbered settings
- Lightweight and compact
- Minimal noise and vibration levels
- Tough durable casing
- On/off switch
- Power on indicator
- Inflation wait indicator
- Audible and visual power fail alarm
- Audible and visual low pressure alarm
- Integrated battery back up for alarms
- Fewer recess areas where dirt can collect
- Hanging hooks
- Snap-fit tube connectors for secure connection
- Easy to use guide on pump



Easy-to-use settings



Seat cushion



Fully alternating modular cells



Durable PU coated cover

Product information

	MATTRESS OVERLAY	MATTRESS REPLACEMENT	SEAT CUSHION	PUMP (CHOOSE REGION)
Order code	ATD002DAR	ATD001DAR	ATD003DAR	UK: ALT01 INDIA: ALT509 KSA: ALT12 USA: ALT03
Dimensions	2040 x 860 x 114 mm (80¾ x 34¾ x 4½")	2040 x 860 x 210 mm (80¾ x 34¾ x 8¾")	455 x 470 x 50 mm (18 x 18½ x 2")	248 x 160 x 116 mm (9¾ x 6¾ x 4½")
Weight				2.7 kg (6 lbs)
Cycle				10 minutes
Power				220-224V, 50 Hz, 14 VA
Electrical safety				IEC60601-1-UK/Euro, IEC60601-1 - ROW
Maximum patient weight	120 kg (265 lbs)	120 kg (265 lbs)	90 kg (198 lbs)	-
Cell material	Polyurethane	Polyurethane	Polyurethane	-
Cover material	Polyurethane coated polyester	Polyurethane coated polyester	Polyurethane coated polyester	-

February 2018 Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. ® and ™ are trademarks belonging to the Arjo group of companies.
© Arjo, 2018

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injury and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Arjo AB · Hans Michelsengatan 10 · 211 20 Malmö · Sweden · +46 10 335 4500

www.arjo.com