



Etac Clean 24"

This variant of Etac Clean is suitable for those who can operate the chair independently and want to maximise their mobility.

Easy to manoeuvre, for the independent user

The Etac Clean 24" wheels are positioned to make it easy to manoeuvre the chair without compromising user safety or comfort. It is easy to transfer to and from the chair thanks to the large, unrestricted seat in front of the wheels. The brake handles are also positioned below the seat so that they do not get in the way during transfer, while remaining easy to reach.

Etac Clean 24" has the same basic benefits as the other members of the Clean family – but allows the user to operate the chair independently.



Accessories

See our wide range of accessories "Accessories for all Clean models" on page 24–25.

Etac Clean 24"

- 55 cm (21¾) with 24" rear wheels, white
- 55 cm (21¾) with 24" rear wheels, lagoon green

Item no.

- 80229276
- 80229216

Materials

- Rear wheel adapter: aluminum
- Rear wheels: polyamide, solid tires
- Chassis: Powder-painted steel tubing
- Seat/back support: Polypropylene
- Arm supports: Polypropylene
- Foot support: Polypropylene
- Castor brackets: Polyamide
- Castors, ball-bearing: Polyamide/TPE

Maintenance



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F if so required.





Self propelled

Easy to operate for the independent user.

Safe and simple transfer

The flat seat and the brake position facilitate lateral transfer.



Accessories



Soft back support

80209225
Adjustable.
Washing at max 60° C (140° F).



Soft pad

Grey: 80209260
Heat-reflecting with a non-slip surface.
Recess width: 17 cm (6¾").
Materials: polyethylene.
Cleaning max max. 60° C (140° F) or with alcohol wipe.




Side Support

80209509
Adjustable in width. Materials: Polyester fabric, Velcro, polyethylene.
Washing at max. 60° C (140° F).
Can not be used in combination with: 80209225.




Soft back pad

Grey: 80209261
Heat-reflecting with a non-slip surface.
Materials: polyethylene. Cleaning max. 60° C (140° F) or with alcohol wipe.




Crossbar

80209246
Attached to the arm supports.
Materials: Aluminium, cellular polyethylene, polypropylene.
Can not be used in combination with: 80209228, 80209237 and 80209226.



Seat cover

80209247
When a complete seat is required.
Materials: Polyurethane.
Cleaning max. 85° C (185° F) or with alcohol wipe.




Chest harness strap

80209442
Adjustable in length. Attach to soft back support (80209225). Material: polyester.



Comfort Kit

80209237
Includes both Comfort Seat Soft and Comfort Cover. Cleaning max. 60° C (140° F) or with alcohol wipe.
Can not be used in combination with: 80209225 and 80209246.




Positioning belt

80209440
Adjustable in length. Attachment with velcro straps.
Material: polyester.



Comfort Cover

80209228
A soft back support with integrated arm support pads. Material: polyethylene. Cleaning max. 60° C (140° F) or with alcohol wipe.
Can not be used in combination with: 80209225 and 80209246.




Safety strap

80209034
Attachments on the rear side of the back support. Always use together with positioning belt. Material: polyester. Can not be used in combination with: 80209225.




Comfort Arm Pads

80209226
Soft and warm, leaving the handles free to ensure a firm grip when standing up.
Material: polyethylene. Cleaning max. 60° C (140° F) or with alcohol wipe. Can not be used in combination with: 80209246.



Cushion, protect

80209245
Air-inflated ROHO cushion with excellent pressure distribution properties. Adjustable air pressure. Recess width: 21 cm (8¼").
Material: neoprene. Cleaning max: 60° C (140° F).



Comfort Seat Soft

Small opening width 15 cm (6"): 80209301
narrow opening width 18 cm (7"): 80209227
A soft contoured seat. Material: polyurethane.
Cleaning max. 85° C (185° F) or with alcohol wipe.

Accessories



Soap dish (does not fit Clean 24")

Grey: 80209266
Attaches to the side of the seat. Built-in shower handle holder.
Material: Polypropylene.
(Does not fit Clean 24")



Bucket holder

Grey: 80209505
Fits Bucket with handle and lid and other buckets with a thick rim. Material: Polypropylene.



Pan

Grey: 80209254
Material: polypropylene.
Pan incl. lid with handle
Grey: 80209255
Material: Polypropylene.



Attachment for low mounting of pan/bucket

Grey: 80209259
Fitted to the pan/bucket holder. Lowers the pan/bucket holder with 6 cm (2½").
Material: Polypropylene.



Pan holder

Grey: 80209257
Fits both pan, bedpan and regular hospital pans with a thin rim.



Disposable pans

80209282
Strong and firm. 50 pcs.



Bucket with handle and lid

80209265
Material: Polypropylene.



Disposable pan holder kit

Grey: 80209229
Fits Disposable pan.

- Clean Height Adjustable only



Flip up foot supports

80209250
Removable, adjustable height. Attachments included.
Materials: powder coated aluminum and polypropylene foot supports.



Elevating leg support


Right: 80209252
Left: 80209253
Attachments included.
Materials: ABS plastic, stainless.

- Clean 49 and 55 cm only




Rear wheel adapter

White: 80209305
Materials: Aluminium, polyamide.



Heightening kit with wheels

White, 4 pcs: 80209511
Lagoon green, 4 pcs: 80209512
Raise the seat by 5 cm (2").
Materials: Polyamide, TPE, aluminium.



75 mm (3") castors

White, front, pair: 80209515
White, back pair: 80209516
Lagoon green, front, pair: 80209513
Lagoon green, back, pair: 80209514
Reduce seat height by 5 cm (2"). Mounted on Clean 44 cm as standard.
Materials: polyamide, TPE, aluminium.