

# Shoprider

## 778SL

ASK FOR THE BRAND **SHOPRIDER**



Active Mobility

RELIABLE MOBILITY GUARANTEED

# Shoprider

## 778SL

**SHOPRIDER®**  
www.shoprider.com.au  
mobility

The SHOPRIDER™ 778SL is a versatile performer from the impressive SHOPRIDER™ range of scooters. The 778SL is loaded with great features to provide the user with outstanding functionality and performance including a powerful motor, large battery capacity, suspension, rotating / sliding seat, modesty panel with dual driving lights & indicators. The reliability of this scooter is further enhanced by the rear wheel drive and electromagnetic brakes.

### SPECIFICATIONS

### 778SL

Length	1230mm
Width	590mm
Height	1250mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	1.2Hp
Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 30km
Wheel Size	260/330mm
Turning Radius	1250mm
Weight (w/o batteries)	76kg

### FEATURES

High Back Seat  
Adjustable Steering Tiller  
Adjustable Width Arm Rests  
Modesty Panel  
Fold Arm Rests  
Sliding Rotating Seats  
Front Shopping Basket  
Variable Seat height  
Battery Level Warning Light  
Front & Rear Suspension  
Front Lights  
Anti Tips Wheels  
**Comprehensive Warranty**

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.


**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING:** Do not exceed 10kph on public roads or paths

#Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

### 778SL



\*(Single user)

Unit 8, 110-120 Silverwater Rd  
SILVERWATER NSW 2128  
Tel: 02 9649 2111  
E: Sales@activemobility.com.au

www.activemobility.com.au



Active Mobility

Updated November 2014