



Etac Clean Height Adjustable

In addition to all the features of the basic model, this variant allows you to adjust the height of the seat.

Adjust the height for a perfect fit

The possibility of adjusting the seat height means that Etac Clean Height Adjustable can be tailor-made to suit the user's needs – for example to make transfers easier or to fit perfectly over the toilet. Adjustment is simple and no tools are required.

Etac Clean Height Adjustable has the same basic benefits as the other members of the Clean family – but is also height-adjustable up to 60 cm, and has a greater range of foot supports and calf supports.

Accessories

See our wide range of accessories "Accessories for all Clean models" on page 24–25.



Etac Clean Height Adjustable

4 lockable wheels, white

Item no.

80229208

Weight

14.2 kg

Maximum user weight

130 kg

Materials

Chassis: Powder coated steel tubing
Seat/Back support: Polypropylene
Arm supports: Polypropylene
Foot support: Polypropylene
Castors with ball, bearings: Polyamide/TPE/Stainless steel tubing
Castor height adjustment tubes: Stainless steel



Maintenance



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F if so required.



Accessories



Soft back support

80209225

Adjustable.

Washing at max 60° C (140° F).



Side Support

80209509

Adjustable in width. Materials: Polyester fabric, Velcro, polyethylene.

Washing at max. 60° C (140° F).

Can not be used in combination with: 80209225.



Crossbar

80209246

Attached to the arm supports.

Materials: Aluminium, cellular polyethylene, polypropylene.


Can not be used in combination with: 80209228, 80209237 and 80209226.



Chest harness strap

80209442

Adjustable in length. Attach to soft back support (80209225). Material: polyester.




Positioning belt

80209440

Adjustable in length. Attachment with velcro straps.

Material: polyester.



Safety strap

80209034

Attachments on the rear side of the back support. Always use together with positioning belt. Material: polyester. Can not be used in combination with: 80209225.



Cushion, protect

80209245

Air-inflated ROHO cushion with excellent pressure distribution properties. Adjustable air pressure. Recess width: 21 cm (8¼"). Material: neoprene. Cleaning max: 60° C (140° F).



Soft pad


Grey: 80209260

Heat-reflecting with a non-slip surface.

Recess width: 17 cm (6¾").

Materials: polyethylene.

Cleaning max max. 60° C (140° F) or with alcohol wipe.



Soft back pad

Grey: 80209261

Heat-reflecting with a non-slip surface.

Materials: polyethylene. Cleaning max. 60° C (140° F) or with alcohol wipe.



Seat cover

80209247

When a complete seat is required.

Materials: Polyurethane.

Cleaning max. 85° C (185° F) or with alcohol wipe.



Comfort Kit

80209237

Includes both Comfort Seat Soft and Comfort Cover. Cleaning max. 60° C (140° F) or with alcohol wipe.

Can not be used in combination with: 80209225 and 80209246.



Comfort Cover

80209228

A soft back support with integrated arm support pads. Material: polyethylene. Cleaning max. 60° C (140° F) or with alcohol wipe.

Can not be used in combination with: 80209225 and 80209246.




Comfort Arm Pads

80209226

Soft and warm, leaving the handles free to ensure a firm grip when standing up.

Material: polyethylene. Cleaning max. 60° C (140° F) or with alcohol wipe. Can not be used in combination with: 80209246.



Comfort Seat Soft

Small opening width 15 cm (6"): 80209301
narrow opening width 18 cm (7"): 80209227

A soft contoured seat. Material: polyurethane.

Cleaning max. 85° C (185° F) or with alcohol wipe.

Accessories



Soap dish (does not fit Clean 24")

Grey: 80209266
Attaches to the side of the seat. Built-in shower handle holder.
Material: Polypropylene.
(Does not fit Clean 24")



Bucket holder

Grey: 80209505
Fits Bucket with handle and lid and other buckets with a thick rim. Material: Polypropylene.



Pan

Grey: 80209254
Material: polypropylene.
Pan incl. lid with handle
Grey: 80209255
Material: Polypropylene.



Attachment for low mounting of pan/bucket

Grey: 80209259
Fitted to the pan/bucket holder. Lowers the pan/bucket holder with 6 cm (2½").
Material: Polypropylene.



Pan holder

Grey: 80209257
Fits both pan, bedpan and regular hospital pans with a thin rim.



Disposable pans

80209282
Strong and firm. 50 pcs.



Bucket with handle and lid

80209265
Material: Polypropylene.



Disposable pan holder kit

Grey: 80209229
Fits Disposable pan.

– Clean Height Adjustable only



Flip up foot supports

80209250
Removable, adjustable height. Attachments included.
Materials: powder coated aluminum and polypropylene foot supports.



Elevating leg support

Right: 80209252
Left: 80209253
Attachments included.
Materials: ABS plastic, stainless.

– Clean 49 and 55 cm only




Rear wheel adapter

White: 80209305
Materials: Aluminium, polyamide.



Heightening kit with wheels

White, 4 pcs: 80209511
Lagoon green, 4 pcs: 80209512
Raise the seat by 5 cm (2").
Materials: Polyamide, TPE, aluminium.



75 mm (3") castors

White, front, pair: 80209515
White, back pair: 80209516
Lagoon green, front, pair: 80209513
Lagoon green, back, pair: 80209514
Reduce seat height by 5 cm (2"). Mounted on Clean 44 cm as standard.
Materials: polyamide, TPE, aluminium.