

# Shoprider

## Rocky 4

ASK FOR THE BRAND SHOPRIDER

NOW WITH  
4 POLE  
MOTOR



ROCKY 4

**SHOPRIDER®**

[www.shoprider.com.au](http://www.shoprider.com.au)

*mobility*

Available Colours



RELIABLE MOBILITY GUARANTEED

# Shoprider

## Rocky 4

**SHOPRIDER**  
www.shoprider.com.au mobility

The SHOPRIDER™ **Rocky 4** is an outstanding model from the extensive range of mobility scooters available from SHOPRIDER. The **Rocky 4** must be seen to be believed, it redefines the concept of reliable mobility. The **Rocky 4** comes with 330mm wheels, Airflow Comfort Seat and all the features you expect from the SHOPRIDER range, fully rotating lockable seat, suspension, adjustable steering tiller and multi position lift up armrests.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

### ROCKY 4

Length	1430mm
Width	680mm
Height	1300mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	10Kph
Battery capacity min	2 x 60 Amp
Maximum range#	up to 50km
Wheel Size	330mm
Turning Radius	1800 mm
Weight (with battery)	146kg

### FEATURES:


Powerful 4 Pole Motor  
Airflow Comfort Seat  
Loop Steering Tiller  
Gas Adjustable Steering Tiller  
Adjustable width arm rests  
Sliding Rotating Seat  
Variable Seat Height  
Battery Power indicator  
Double Head Lights  
Rear Lights & Indicators  
Hazard Lights  
Front & Rear Bump Pads  
High & Low Speed  
Shopping Basket  
**Comprehensive Warranty**

**WARNING:** Do not exceed 10kph on public roads or paths

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**SHOPRIDER**  
www.shoprider.com.au mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at [www.shoprider.com.au](http://www.shoprider.com.au)