

DRY FLOATATION® Mattress Overlay



Non-powered for
easy portability

Drainage holes in base
to draw away fluids

Flame resistant, durable
neoprene construction
for easy cleaning

The ROHO DRY FLOATATION Mattress Overlay is an adjustable, zoned, air filled pressure reducing surface suitable for patients at all levels of risk.

- 4 independently adjustable, interconnected zones
- flame resistant, durable neoprene construction for easy cleaning
- non-powered for easy portability
- low friction and shear environment
- midweight overlay (3.6 kg per section; 14.4 kg complete system)
- drainage holes in base to draw away fluids
- suitable for users of any weight



ROHO DRY FLOATATION products are made of soft, flexible interconnected air cells designed to conform to each user's individual body shape and self-adjust with every movement, enhancing blood flow and supporting posture.

Individual, interconnected cells move independently, allowing each cell to twist, turn, bend and adapt precisely to the contours and anatomy of the user providing a total contact environment. Cells mimic the properties of water by evenly distributing pressure across the cushioning surface and away from peak pressure areas that are most at risk for developing pressure injuries. The low surface tension allows for immersion without tissue deformation, minimising the risk of skin breakdown.



The DRY FLOATATION Mattress Overlay system is comprised of four interconnected zones (each 510 x 860 mm) made from independent, neoprene rubber air cells that encourage air movement and ventilation.

DRY FLOATATION ® Mattress Overlay

Enclosure cover
included



4 independently adjustable,
interconnected zones



The 80 mm high cells can be adjusted per segment to individual patient body shapes and contours, with no patient weight restrictions.



Providing both friction and shear relief, the durable neoprene rubber is free of latex, flame-resistant and easy to clean. Drainage holes in the overlay base help to pull fluids away from incontinent patients.

Includes a free ROHO DRY FLOATATION Enclosure Cover* - a fluid resistant, fully enclosed mattress overlay cover for increased infection control and incontinence protection. Straps secure the cover to a standard single bed mattress, while carry handles help with product transport. Additional covers are also available as an optional accessory.

Also comes complete with a hand inflation pump, repair kit and operating instructions.

| ITEM | Mattress Overlay | Single Section | Enclosure Cover |
|-------------|------------------|------------------|------------------|
| CODE | R-MATT SYSC | R-MATT RESS-SECT | R-COV-MR2PCSTD-N |
| LENGTH | 1905 mm | 510 mm | |
| WIDTH | 860 mm | 860 mm | |
| HEIGHT | 85 mm | 85 mm | |
| CAPACITY | Unlimited | Unlimited | |
| WEIGHT | 14.4 kg | 3.6 kg | |
| CONFORMS TO | ARTG 111257 | | |
| WARRANTY | 2 year | | |

* Enclosure cover for the Mattress Overlay only (not available for the single section)

