



## Select Series

The QUADTRO SELECT® and CONTOUR SELECT® cushions feature four independent compartments of air cells for an individualized fit, each controlled by our exclusive ISO-FLO® Memory Control.

ROHO's ISO-FLO® Memory Control, located at the front of all Select Series cushions, allows air to flow through the cushion in two ways:



**ISO-FLO Memory Control shown in the OPEN position.**

Air can flow freely throughout the entire cushion, similar to a ROHO Single Compartment Cushion.



**ISO-FLO Memory Control shown in the CLOSED position.**

Air becomes locked into four separate quadrants, offering front-to-back and side-to-side positioning capabilities.

Every Select Series cushion sets the standard for skin and soft tissue protection, without compromising stability, positioning, comfort or ease of use.

Featuring DRY FLOATATION® air cells, the Select Series assists in pressure injury prevention and treatment.

## QUADTRO SELECT® CUSHIONS



### HIGH PROFILE®

The tall interconnected air cells envelop your body, delivering maximum immersion for the ultimate in skin protection, positioning, stability and comfort.

#### IDEAL FOR PEOPLE WHO

Require accommodation or correction of front-to-back or side-to-side pelvic asymmetries, up to approx. 3 in. (7.5 cm)

Require pelvic and lower extremity positioning

#### DIMENSIONS

Width: 12.75 in. to 26.25 in. / 32.5 cm to 67 cm  
Depth: 11.75 in. to 67 in. / 30 cm to 51.5 cm  
Height: 4.25 in. / 10.5 cm

#### LIMITED WARRANTY

Cushion: 36 months  
Cover: 6 months

#### WEIGHT LIMIT

No weight limit

#### MEDICARE CODES

E2624, E2625



### MID PROFILE™

The mid-sized interconnected air cells deliver the ideal combination of customizable positioning and stability.

#### IDEAL FOR PEOPLE WHO

Require accommodation or correction of front-to-back or side-to-side pelvic asymmetries, up to approx. 2 in. (5 cm)

Require pelvic and lower extremity positioning

#### DIMENSIONS

Width: 14.5 in. to 19.5 in. / 36.5 cm to 49.5 cm  
Depth: 13.5 in. to 18.5 in. / 34.5 cm to 47 cm  
Height: 3.25 in. / 8.5 cm

#### LIMITED WARRANTY

Cushion: 36 months  
Cover: 6 months

#### WEIGHT LIMIT

No weight limit

#### MEDICARE CODES

E2624, E2625



### LOW PROFILE®

The interconnected air cells deliver sitting stability. It could be used as a secondary cushion, on an office chair, in adaptive sports equipment, in a car or anywhere else you may need cushioning and protection.

#### DIMENSIONS

Width: 11 in. to 24.25 in. / 28 cm to 62.5 cm  
Depth: 11.75 in. to 22 in. / 30 cm to 56 cm  
Height: 2.25 in. / 5.5 cm

#### LIMITED WARRANTY

Cushion: 36 months  
Cover: 6 months

#### WEIGHT LIMIT

No weight limit

#### MEDICARE CODES

E2624, E2625

## CONTOUR SELECT®



Designed for positioning and posture, this cushion helps to keep you centered in the cushion in an upright pelvic posture.

- Provides leg alignment and pelvis stabilization
- Provides further individualization of side-to-side and front-to-back positioning with ISOFLO® Memory Control
- Stability
- Positioning

#### DIMENSIONS

Width: 14.5 in. to 24.25 in. / 36.5 cm to 62.5 cm  
Depth: 15.25 in. to 20.25 in. / 38.5 cm to 51.5 cm  
Height: Approx. 2.25 in. to 4.25 in. / 5.5 cm to 10.5 cm  
Range due to varying heights of cells

#### LIMITED WARRANTY

Cushion: 36 months  
Cover: 6 months

#### WEIGHT LIMIT

No weight limit

#### MEDICARE CODES

E2624, E2625